

Waste Separator

PHS Wastemanagement offers a complete service for the disposal of waste from Dürr and other leading brand amalgam separator systems, to ensure your compliance with waste management regulations.

Product Information

Amalgam Separator cup **Small, Medium, Large**

EWC Code **18 01 10**

We also offer a disposal service for non-Dürr systems, please ask your PHS representative.

Instructions for use:

Users should follow manufacturer's instructions and ensure that they have the correct separator cup for their system.



Amalgam Sludge Disposal

The PHS Wastemanagement Amalgam Sludge container is a convenient and robust unit designed for the secure storage of amalgam sludge. To help protect staff and patients from the potentially harmful effects of mercury inhalation, the container has an integral mercury suppressant foam.

Product Information

Capacity	6L
Size	260mm x 200mm (dia)
EWC Code	18 01 10

Instructions for use:

Simply store used amalgam sludge in this container until ready for collection by your PHS representative.

